



Easy Setup for Radiology Centres

Regulatory Cleared Service Delivery – FDA, CE Mark, and TGA

Following an easy setup process, Resonance Health provides a range of clinical image analysis services of MRI scans via a unique quality assured service delivery model:

- ✓ FerriScan® - the global gold standard for the measurement of liver iron concentration (LIC)
- ✓ Cardiac T2* - assessment of cardiac iron (FerriScan add-on)
- ✓ HepaFat-Scan® - measurement of volumetric liver fat fraction (VLFF)

Image data is securely transferred to Resonance Health's central ISO certified Service Centre, where it is analysed by our team of experts. Result reports are then made available for download by authorised Radiology staff.



- ✓ Services are charged per analysis only - there are no subscription or licence fees
- ✓ Radiology Centres are not required to purchase any additional hardware or software
- ✓ The imaging protocols require no contrast agents and can be established on most 1.5 Tesla scanners

Easy Setup Process



1. **Paperwork:** Scanner details collected, scanner compatibility confirmed, and contractual arrangements established - *there is no minimum use requirement*



2. **Secure Image Transfer:** Account created for secure, HIPAA compliant web-based FAST system - *for image data transfer, analysis tracking, and result downloading*



3. **Quality Control and Standardisation Checks:** Radiology Centre scans FerriScan Phantom Pack (FerriScan only) or test subject (Cardiac T2* and HepaFat-Scan) - *data submitted via FAST to Resonance Health for quality control checks*



4. **Patient Scanning Confirmation:** FerriScan Phantom Pack or test subject image data analysed at Resonance Health central Service Centre - *confirmation that patient scanning can commence or technical support to ensure quality of results*

To arrange setup at your Centre please contact info@resonancehealth.com